J E J O U E

Ami

Introducing Ami.

The pleasurable way to tone your PC muscle.

Ami is unique.

Its innovative three-step system allows you to progress at your own pace.

Ami is pleasure.

As the hidden weights move, teasing as you tone, exercise becomes enjoyment.

Ami is a friend.

Helping you experience greater control and more intense orgasms.

Why use Ami?

Pleasure for you

Using Ami regularly will tone and strengthen your PC muscle, so you can enjoy enhanced sensitivity and more intense, easy-to-achieve orgasms. What's more, each Ami contains little balls that rumble around inside you, creating a discreet, arousing sensation even while you're working out.

Pleasure for him

If you have a male partner, Ami can lead to a whole new sexual experience when you use your improved PC muscle to grip and tease his penis during intercourse.

Health

We might not like to talk about it, but many women experience a degree of urinary incontinence at some point in their lives. With Ami, you can strengthen the PC muscle to improve bladder control.

Pregnancy

Using Ami before you fall pregnant can protect your vagina by preparing the PC muscle for delivery. After birth, Ami can speed up vaginal recovery, and the enhanced tightness you'll achieve will increase sensitivity for you and your partner during sex.

Menopause

Strengthening your PC muscle can enhance your natural vaginal lubrication, sexual response and urinary control during the physiological changes of menopause.

Your PC muscle

The pubococcygeus (PC muscle) is the hammock-like muscle that stretches from your pubic bone to coccyx, supporting your pelvic organs and internal genitalia. A strong PC muscle is the key to more intense orgasms and greater urinary control, as well as helping during childbirth.

Kegel exercises

Kegel exercises (named after the American gynecologist, Dr. Kegel, who popularised them in the 1970s), are designed to tone and strengthen your PC muscle. With Ami, these exercises become simple, easy and pleasurable.

The three Amis

Ami isn't just one Kegel exerciser, but three – making it a totally unique way to tone up your pelvic floor. We've taken cues from the best medical and pleasure products on the market to create a personal workout program that's a whole lot more enjoyable than the gym.

Ami 1

A large, soft and lightweight single ball. If you're new to Kegel exercise products or you're out of practice, we recommend starting here.

Ami 2

A medium-sized double ball that's still soft yet slightly heavier. Try Ami 2 if exercising with Ami 1 doesn't challenge you enough.

Ami 3

A small, hard and heavy double ball – this is the most advanced Ami for experienced users.

Using Ami

Before enjoying Ami, you'll need to locate your PC muscle. You can easily find it by trying to stop the flow of urine while going to the loo. The muscle you're contracting is your PC muscle.

Now you've found it, choose your Ami. Beginners should start with Ami 1, and you'll gradually be able to progress to Ami 2 and Ami 3 as your muscle strength improves. Coat your Ami with generous amounts of water-based lubricant before you begin exercising.

There are many different Kegel exercises, and some may suit you more than others. Try the popular ones below or experiment with other techniques, such as wearing Ami around the house to enjoy the pleasant vibrating sensation it creates.

Exercise 1: Something simple

Relax down onto your back on a bed or the floor and insert your Ami. Squeeze your PC muscle for five seconds to grip Ami, then release for five seconds. Repeat this 10 times.

Take a one-minute break, then repeat this again. Next, take a two-minute break before doing a third and final set of 10.

Perform this daily, and when it gets too easy, try three sets of 15 instead of 10.

Exercise 2: A little more advanced

Sit on the edge of a bed or chair and insert the Ami of your choice. Squeeze your PC muscle, then use the string to gently pull your Ami out for five seconds, then relax for five seconds. Repeat this 10 times.

Give yourself a one-minute break before repeating this again. Next, take a two-minute break before your last set of 10.

Do this every day, and when you feel it's too simple, increase the number of repetitions to three sets of 15 instead of 10.

Spec:

Materials:

100% body safe silicon and Polycarbonate plastic

Weight:

Ami 1 – 47 g / 1.65 oz

Ami 2 – 78 g / 2.75 oz

Ami 3 – 106 g / 3.73 oz

Length inc. cord:

Ami 1 – 111 mm / 4.37 inches

Ami 2 – 156 mm / 6.14 inches

Ami 3 – 144 mm / 5.66 inches

Length excluding cord:

Ami 1 – 36 mm / 1.41 inches

Ami 2 – 80 mm / 3.14 inches

Ami 3 – 70 mm / 2.75 inches

Maximum diameter:

Ami 1 – 36 mm / 1.41 inches

Ami 2 – 32 mm / 1.25 inches

Ami 3 – 28 mm / 1.10 inches

Lubrication:

We recommend using only water-based lubricants with Ami

For cleaning:

We recommend washing Ami in warm water and antiseptic soap. However if you would like to boil Ami, we recommend that you do so only for 1-2 minutes maximum, in order to prolong the usable life of Ami.

Note:

Ami is not a medical device. Neither these statements nor this product have been examined, approved or evaluated by the FDA. This product is not intended to diagnose, treat, cure or prevent any disease.

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IMPORTANT SAFETY INFORMATION
READ ALL INSTRUCTIONS BEFORE USE.
DO NOT DISCARD THESE INSTRUCTIONS.

To reduce the risk of injury and ensure your enjoyment of AMI:

1. Use AMI only in accordance with these instructions.
2. AMI is not a child's toy. Please keep away from children.
3. DO NOT operate AMI if the outer layer of silicone has been damaged.
4. AMI serves exclusively for body massage and intimate stimulation of erogenous zones. DO NOT use AMI for other uses.
5. It is important to note that pathogens could be transmitted between people by sharing your AMI. To protect against the risk of transmitting pathogens, do not share your AMI with other people.
6. DO NOT use AMI if you encounter any discomfort during or after use.
7. DO NOT use AMI as a substitute for medical attention.
8. DO NOT use AMI on sensitive skin areas or in the presence of poor circulation.
9. AMI is not intended for use by persons with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of AMI by a person responsible for their safety.
10. Always clean your AMI after use to reduce the risk of infections.
11. AMI is not intended as a medical appliance.

If you have concerns regarding your health, consult your doctor before using AMI.

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